

<i>NO</i>	<i>Test name</i>	<i>Sample type</i>	<i>S.Vol (ml)</i>	<i>Assay method</i>	<i>Storage</i>	<i>Transport</i>	<i>Preparation before sampling</i>	<i>Reject</i>	<i>Bilirubin</i>	<i>Hemolysis</i>	<i>Lipemia</i>	<i>Storage conditions after test performance</i>
1	T3	S	0.5	CLIA	8 hours at 2-8°C, 48 hours at -20°C	at 2-8°C, - 20°C	Notice to hormonal treatment & gestation	Hemolysis, Lipemic, ictric	No effect Up to 200 mg/l bilirubin	No effect Up to 250 mg/dl hemoglobin	No effect Up to 1000 mg/dl Triglycerides	1 months at -20
2	T4	S	0.5	CLIA	8 hours at 2-8°C, 48 hours at -20°C	at 2-8°C, - 20°C	Notice to hormonal treatment & gestation	Hemolysis, Lipemic, ictric	No effect Up to 200 mg/l bilirubin	No effect Up to 500 mg/dl hemoglobin	No effect Up to 1000 mg/dl Triglycerides	1 months at -20
3	TSH	S,PE,PH	0.5	CLIA	48 hours at 2- 8°C, 30 days at - 20°C	at 2-8°C, - 20°C	Notice to hormonal & chemotherapy Drugs	Hemolysis, Lipemic, ictric	No effect Up to 400 mg/l bilirubin	No effect Up to 500 mg/dl hemoglobin	No effect Up to 1000 mg/dl Triglycerides	1 months at -20
4	Amilase	S,PE,PH,urine	0.5	Enzymatic&colorimetric	Serum&Plasma:7 Days at 2-8°C, 30 days at -20°	at 2-8°C, - 20°C	8 hrs Fasting recommended	Hemolysis, Lipemic, ictric	No effect Up to 30 mg/dl bilirubin	effect Up even in small amounts	No effect Up to 800 mg/dl Triglycerides	1 months at -20
5	Lipase	S,PH,urine	0.5	colorimetric	4 Days at 2-8°C, 6 months at -20°	at 2-8°C, - 20°C	Fasting recommended	Hemolysis, Lipemic, ictric	No effect Up to 300 mg/dl bilirubin	No effect Up to 200 mg/dl hemoglobin	No effect Up to 800 mg/dl Triglycerides	1 months at -20
6	Gamma G.T	S,PE	0.5	colorimetric	30 Days at 2-8°C, 6 months at - 20°,7 days Ambient	at 2-8°C, - 20°C	2 hrs Fasting recommended	Hemolysis, Lipemic, ictric	No effect Up to 40 mg/dl bilirubin	No effect Up to 100 mg/dl hemoglobin	No effect Up to 800 mg/dl Triglycerides	1 months at -20

WB : Whole blood      S : Serum      U : Urine      BF : Body fluids      H : Heparin

D: Defibrinate      PH: Heparinized plasma      PE: Plasma (EDTA)      PC: Plasma (Citate)